

DISPOSABLE ISOLATION GOWN

All gowns are produced with the best commercial efforts to adhere to standards set in place by the FDA, CDC and OSHA, though they have been slightly modified to address specific needs of the COVID-19 response effort.

STYLE GUIDE

Wrist styles

Thumb hook



Knit cuff

Neck style

Over-the-head



Front / back styles

Open back




DISPOSABLE ISOLATION GOWN

MATERIAL

Spunbond Polypropylene, non-woven single layer fabric

<u>Properties</u>	<u>Units</u>	<u>Target</u>	<u>Max</u>	<u>Min</u>
Basis Weight	Oz/ SY	1.5	1.66	1.34
MD Grab Tensile	LBS/IN	31.5	41.9	21.1
MD Grab Elongation	%	57	81	33
CD Grab Tensile	LBS/IN	25.4	33.8	17
CD Grab Elongation	%	75	113	44
MD Trap Tear	LBS	17.3	24.5	10.1
CD Trap Tear	LBS	12.6	18.6	6.6
Air Permeability	M3/M2/MIN	114	165	63
Water Flow Rate	GAL/MIN/FT ²	15.2		



DISCLAIMER, BUYER'S RISK, PLEASE READ

The product noted above has been produced by Artisan's Asylum using commercially reasonable efforts to source materials and provide finished product that complies with standards acceptable to the Buyer. The Buyer acknowledges that they have reviewed the product offering from Artisan's Asylum and that no guarantees are provided to the Buyer for the specific quality, size, consistency, material, availability, delivery, or efficacy of the product. The Buyer is purchasing the product at their own risk. This product may contain latex.

Artisan's Asylum will make commercially reasonable efforts to fulfill the order. Due to the COVID-19 crisis, material supplies, and crisis situation globally, payment must be made to Artisan's Asylum in advance and there are no refunds for purchases made. The Buyer accepts all potential of loss.